

Screen Quest™ Colorimetric ELISA cAMP Assay Kit

 Catalog number: 36370, 36371
 Unit size: 1 plate, 10 plates

Component	Storage	Amount (Cat No. 36370)	Amount (Cat No. 36371)
Component A: cAMP Standard	Freeze (< -15 °C), Minimize light exposure	1 vial (33 µg)	1 vial (33 µg)
Component B: Assay Buffer	Refrigerated (2-8 °C)	1 bottle (20 mL)	1 bottle (100 mL)
Component C: HRP-cAMP Conjugate	Freeze (< -15 °C), Minimize light exposure	1 vial	5 vials
Component D: 10X Wash Solution	Refrigerated (2-8 °C)	1 bottle (10 mL)	1 bottle (100 mL)
Component E: Cell Lysis Buffer	Refrigerated (2-8 °C)	1 bottle (10 mL)	1 bottle (100 mL)
Component F: Anti-cAMP Ab Coated 96-Well Plate	Refrigerated (2-8 °C), Minimize light exposure	1 plate	10 plates
Component G: Amplitude® Green	Refrigerated (2-8 °C), Minimize light exposure	1 bottle (10 mL)	1 bottle (100 mL)

OVERVIEW

Screen Quest™ Colorimetric ELISA cAMP Assay Kit is designed for measuring adenosine 3',5'-cyclic monophosphate (cAMP), an important second messenger in intracellular signal transduction. Monitoring cAMP levels is commonly used for screening GPCR agonists and antagonists, especially for receptors coupled to adenylate cyclase-regulated signaling pathways. The assay is based on competition between HRP-labeled cAMP and unlabeled free cAMP from the sample. In the absence of sample cAMP, HRP-cAMP binds to the anti-cAMP antibody. When unlabeled cAMP is present, it competes with HRP-cAMP for antibody binding, reducing the amount of HRP-cAMP bound in the assay. Amplitude® Green is used as the colorimetric HRP substrate to generate the measurable signal, and the assay can be performed in 96-well or 384-well microplate formats. Screen Quest™ Colorimetric ELISA cAMP Assay Kit is suitable for cAMP quantitation, GPCR signaling studies, adenylate cyclase pathway analysis, and biochemical or cell-based screening assays.

AT A GLANCE
Protocol Summary

1. Prepare samples
2. Add 75 µL/well of cAMP standard or test samples into the anti-cAMP coated 96-well plate
3. Incubate at room temperature for 5-10 mins
4. Add 25 µL/well of 1X HRP-cAMP Conjugate
5. Incubate at room temperature for 3 hours
6. Wash 4 times with 200 µL/well Washing Buffer
7. Add 100 µL/well of Amplitude® Green
8. Incubate at room temperature for 1 to 3 hours
9. Monitor absorbance increase at 405, 650 or 740 nm

Important Note

Do not freeze Anti-cAMP Ab Pre-coated 96-well plate (Component F), store it at 4°C. Allow all the kit components to warm to room temperature before using them. Some material might be stick to the vial cap during the shipment. Briefly centrifuge the vial to collect all the content.

KEY PARAMETERS
Absorbance microplate reader

Absorbance 405, 650, or 740 nm
 Recommended plate Clear plate (Component F)

CELL PREPARATION

For guidelines on cell sample preparation, please visit <https://www.aatbio.com/resources/guides/cell-sample-preparation.html>

PREPARATION OF STOCK SOLUTIONS

Unless otherwise noted, all unused stock solutions should be divided into single-use aliquots and stored at -20 °C after preparation. Avoid repeated freeze-thaw cycles

cAMP stock solution (100 µM)

Add 1 mL of Assay Buffer (Component B) to the vial of cAMP Standard (Component A).

HRP-cAMP conjugate stock solution (50X)

Add 55 µL (Cat. # 36370) or 110 µL (Cat. # 36371) of Assay Buffer (Component B) into the vial of HRP-cAMP Conjugate (Component C).

Note: The unused 50X HRP-cAMP conjugate stock solution should be divided into single use aliquots and stored them at -80 °C after preparation and used within 4 weeks.

Washing solution (1X)

Add 1 mL of 10X Wash Solution (Component D) to 9 mL distilled water.

PREPARATION OF STANDARD SOLUTIONS

For convenience, use the Serial Dilution Planner: <https://www.aatbio.com/tools/serial-dilution/36370>

cAMP standard

Make 1:10, 1:100 and 1:3 serial dilutions of cAMP standards in Assay Buffer (Component B) to have 10,000, 100, 30, 10, 3, 1, 0.3, 0.1, 0.03, 0.01, and 0.003 nM cAMP diluted solutions. Store on ice or 4°C.

PREPARATION OF WORKING SOLUTION
HRP-cAMP Conjugate working solution

Make 1:50 dilution with Assay Buffer (Component B) to have 1X HRP-cAMP conjugate working solution before use. Store it on ice or 4°C.

Note: 25 µL of 1X HRP-cAMP conjugate working solution is enough for one assay point; prepare appropriately volume for single use only.

SAMPLE EXPERIMENTAL PROTOCOL
Prepare samples

1. Treat cells as desired: The following is an example of HeLa cells treated with Forskolin to induce cAMP in a 96-well plate format: Aspirate off cell growth medium, add 100 μL /well 100 μM Forskolin in Hanks and 20 mM Hepes buffer (HHBS), incubate in a 5% CO_2 , 37°C incubator for 15 minutes. Aspirate off cell solution after the incubation, add 100 μL /well of Cell Lysis Buffer (Component E), and incubate at room temperature for another 10 minutes. This cell lysate can be assayed directly or after diluted in Assay Buffer (Component B).

Note: Each cell line should be evaluated on an individual basis to determine the optimal cell density. Cells may be seeded the day before or on the day of the experiment depending upon the cell type and/or the effect of the test compounds.

2. Tissue Samples: It is important to rapidly freeze tissues after collection (e.g., using liquid nitrogen) due to quick metabolism of cyclic nucleotides in tissue. Weigh the frozen tissue and add 10 - 20 μL /mg of cell lysis buffer. Homogenize the sample on ice. Spin at top speed for 5 minutes and collect the supernatant. The supernatant may be assayed directly.
3. Urine, Plasma and Culture Medium Samples: Urine and plasma may be tested directly with 1:200 to 1:1000 dilutions in 1X Lysis Buffer. Culture medium can also be tested with 1:10 to 1:200 dilutions in Lysis Buffer.

Note: RPMI medium may contain > 350 fmol/ μL cAMP.

cAMP assay

1. All the assay wells will be prepared in the following orders: A) cAMP standards, control, or test samples; B) HRP-cAMP Conjugate.
2. Add 75 μL /well of the cAMP diluted standard solution and test samples into each well of the anti-cAMP Ab coated 96-well plate (Component F). We recommend duplicating the assays for each standard and testing sample. Incubate at room temperature for 5 to 10 minutes.
3. Add 25 μL /well of 1X HRP-cAMP Conjugate working solution. Incubate at room temperature for 3 hours by placing the plate on shaker.
4. Aspirate plate contents, and wash 4 times with 200 μL /well of 1X wash solution.
5. Add 100 μL /well of Amplite® Green (Component G) into each well, and incubate at room temperature for 60 mins to 3 hours, protected from light.
6. Monitor the absorbance increase at 405 nm, 650 nm, or 740 nm using an absorbance plate reader.

EXAMPLE DATA ANALYSIS AND FIGURES




Figure 1. cAMP dose response was measured with Screen Quest™ Colorimetric ELISA cAMP Assay Kit in a clear 96-well plate with a SpectraMax microplate reader. The Absorbance can be read at 405 nm (blue line), 650 nm (red line) or 740 nm (Green line), the data in figure B are from the incubation with Amplite® Green for 3 hours.

DISCLAIMER

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email info@aatbio.com if you have any questions.