

CBCP6 Antibody [#C15302]

Catalog Number:	C15302
Concentration:	1mg/ml
Swiss-Prot No.:	Q5VU57
Other Names:	ATP/GTP-binding protein-like 4; CBPC6; CCP6; cytosolic carboxypeptidase 6; EC 3.4.17.-
Storage/Stability:	Store at -20°C/1 year
Form of Antibody:	Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human CBCP6.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Specificity:	CBCP6 antibody detects endogenous levels of total CBCP6 protein.
Reactivity:	Human
Applications:	WB: 1:500~1:1000 ELISA: 1:20000
References:	Ota T., Nat. Genet. 36:40-45(2004). Gregory S.G., Nature 441:315-321(2006).