# SFT2B Antibody［#C19499］

| | |
|---|---|
| **Catalog Number:** | C19499 |
| **Concentration:** | 1mg/ml |
| **Swiss-Prot No.:** | O95562 |
| **Other Names:** | dj747l4.c1.2; flj34085; rp4-747l4.1; sft2 domain containing 2; sft2b; sft2d2; unq512 |
| **Storage/Stability:** | Store at -20°C/1 year |
| **Form of Antibody:** | Rabbit IgG in phosphate buffered saline (without $Mg^{2+}$ and $Ca^{2+}$), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| **Immunogen:** | The antiserum was produced against synthesized peptide derived from internal of human SFT2B. |
| **Purification:** | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| **Specificity:** | SFT2B antibody detects endogenous levels of total SFT2B protein. |
| **Reactivity:** | Human (Identities = 100%, Positives = 100%); Mouse (Identities = 100%, Positives = 100%); Rat (Identities = 100%, Positives = 100%) |
| **Applications:** | WB: 1:500~1:1000 ELISA: 1:40000 |
| **References:** | Rhodes S., Submitted (JAN-1999) to the EMBL/GenBank/DDBJ databases.<br>Gregory S.G., Nature 441:315-321(2006).<br>The MGC Project Team, Genome Res. 14:2121-2127(2004). |

**©2011 by AAT Bioquest, Inc.**, 520 Mercury Drive, Sunnyvale, CA 94085. Tel: 408-733-1055
**Ordering:** sales@aatbio.com; 800-990-8053 or 408-733-1055; Fax: 408-733-1304
**Technical Support:** support@aatbio.com; 408-733-1055